


P/N: 25041-C



Your keys to AssayMAP technology are the 24 PA50 cartridges contained in the 2 tubes attached to this card. These cartridges use the same Protein A chemistry employed in conventional affinity HPLC assays, but offer improved productivity by means of a 100-fold increase in throughput while using standard laboratory equipment and your current validated reagents.


 **AssayMAP**  
\_breakthrough



24 cartridges



AssayMAP is technology patented by BioSystem Development. AssayMAP® is a registered trademark of BioSystem Development, Madison, WI. AssayMAP PA cartridges are made with POROS® MabCapture® Media from Life Technologies Corporation. POROS® and MabCapture® are registered trademarks of Life Technologies Corporation.



802 Deming Way, Suite 172  
Madison, WI 53717 USA  
P: 608.836.5552 F: 608.836.5558  
[info@biosystemdevelopment.com](mailto:info@biosystemdevelopment.com)



breakthroughput

P50035 — AssayMAP® PA50 Starter Kit (24 cartridges)